


# ZOLL AED PLUS

OK?



LIFT TO OPEN



ON  
OFF



# What good is an AED that only works half of the time?

## It's a real rescue...

A sudden collapse—the reason you have an AED. You attach electrode pads to the collapsed victim's bare chest and your AED begins analyzing the heart rhythm. Soon, it's going to tell you to shock or not.

## Half of the time, your AED will say "No shock advised"

Many people don't realize that in about half of all rescues performed by responders, No shock is advised\*.

And, when no shock is advised, providing CPR to the victim of cardiac arrest is the most important factor affecting survival. It's even possible that, with high-quality CPR, the next heart rhythm analysis will tell you to push the shock button.

## But will your AED help you do CPR?


An AED from ZOLL® equipped with Real CPR Help® and CPR-D padz attached will. Other AEDs can't see when you're doing CPR and can't help you do proper compressions at the correct depth and rate. An AED that says "No shock advised," but then doesn't help you do CPR, is really only working half the time.

## Only working half the time?

That's right. By not "seeing" your compressions, AEDs without Real CPR Help will not help you provide the high-quality CPR needed for the best chance of survival.

With Real CPR Help, you can be confident that when your AED says "Start CPR," you'll get real-time feedback that helps you with depth and rate of compressions because you'll see what's happening when providing chest compressions.

\* ML Weisfeldt, et. al. for the Resuscitation Outcomes Consortium (ROC) Investigators, Ventricular Tachyarrhythmias after Cardiac Arrest in Public versus at Home, *The New England Journal of Medicine* 2011;364:4.




# ZOLL®

ZOLL Medical Corporation  
269 Mill Road  
Chelmsford, MA 01824

978-421-9655  
800-348-9011  
[www.zoll.com/aedplus](http://www.zoll.com/aedplus)

©2011 ZOLL Medical Corporation. All rights reserved. AED Plus, Real CPR Help and ZOLL are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners. 9656-0161 0807100





ADP PLUS  
BATTERY BOARD  
COPYRIGHT © 2008  
T10Z



DURACELL  
ULTRA  
123 3.0 V  
LITHIUM 3.0 V  
CR17345

ASSEMBLED IN  
USA

X11C618639

